

Summary

In order, the top three consumables of an FTIR are the desiccant, IR source, and HeNe laser.

The HeNe laser is a critical component of the Interferometer and is used to calibrate the instrument as well as control the moving mirror.

Replacement of the HeNe laser is often needed due to:

- High laser gains after alignment (pegged at maximum of 240).
 - Scan issues resulting to no scanning, intermittent bad scans, or noisy data.
 - No visible laser light in sample compartment or from the laser head (dead laser).
 - Old lasers (from the manufacturing date shown on the head) of more than 7 years will often show one or more of these top three symptoms.
-

Removing the old laser

1. Turn off the FTIR.
2. Remove back Panel (single screw).



Removing the old laser (continued)

3. Unplug the laser power connector on the laser module.



4. Using a large Philips screw driver loosen screws, and pull module out.



Installing the new laser

1. Remove tape on laser (protects the lens from dust/handling).
2. Reverse the same sequence as above in the removal process.
3. The laser can only go in one way. On the bottom side of the module there are dowel pins that assure the laser will be seated for proper alignment.
4. Power the FTIR up.
5. Open the Omnic Application.
6. Go into Experiment Setup, then the Diagnostics tab.
7. Hit Align and verify the alignment has completed.
8. When the alignment is complete you should see an interferogram of some amplitude as shown to the right.
9. Exit Omnic and go into Local Disk C
10. Program Files->Omnic-> and run the application called "BenchStatusReport (Running this status report application will report all settings in a TXT file



